

# EasyFlash Columns

Chromservis has introduced Flash chromatography columns for high performance purifications. New **EasyFlash** columns offer high efficiency for a fast separation and isolation of impurities, reaction products and intermediates.

We offer wide range from 4 g up to 1600 g **EasyFlash** columns with various silicagel types – spherical and irregular.




**EASY  
FLASH**

## Flash sorbent types:

Description	Type	50 µm	30 µm	25 µm	15 µm	32/63 µm
EasyFlash Silica IR	Irregular	•				
EasyFlash Silica 68	Spherical	•		•	•	
EasyFlash Silica 50	Spherical	•	•		•	
EasyFlash C18 High Performance	Spherical	•	•		•	
EasyFlash C18-AQ	Spherical		•		•	
EasyFlash Alumina Basic	Spherical					•
EasyFlash Alumina Neutral	Spherical					•
EasyFlash NH2	Spherical	•	•		•	
EasyFlash SAX	Spherical	•				
EasyFlash SCX	Spherical	•				

## Technical parameters of the EasyFlash columns:

Size	CV [mL]	Typical flow rate [mL/min]	Length [mm]	ID [mm]
4 g	5	5	68	12
12 g	19	15	84	21
25 g	32	15	133	21
40 g	48	25	135	27
80 g	102	35	205	31
120 g	153	45	224	36
220 g	269	130	153	60
330 g	405	130	226	60
800 g	1078	220	341	78
1600 g	2170	380	385	104

For more details please see here 

Ask for your quote by email: [inquiries@chromservis.eu](mailto:inquiries@chromservis.eu)

For more details please log on our web: [www.chromservis.eu](http://www.chromservis.eu)

